

## Hei-VAP Core Handlift G1

### Technical datasheet



|                                                  |                                           |
|--------------------------------------------------|-------------------------------------------|
| <b>hei.techdata.Standard Supply Voltage</b>      | 230 V 50/60 Hz or 115 V 50/60 Hz          |
| <b>hei.techdata.Power consumption (W)</b>        | 1400                                      |
| <b>hei.techdata.Device fuses</b>                 | T8 H (230 V) or T15A H (115 V)            |
| <b>hei.techdata.Protection Class (IEC 61140)</b> | 1                                         |
| <b>hei.techdata.Protection Class (IEC 60529)</b> | Device: IP 20 / Cable heating bath: IP 67 |