

## Hei-VAP Expert Technical datasheet



<b>hei.techdata.Standard Supply Voltage</b>	230 V 50/60 Hz or 115 V 50/60 Hz
<b>hei.techdata.Power consumption (W)</b>	1400
<b>hei.techdata.Device fuses</b>	T8 H (230 V) or T15A H (115 V)
<b>hei.techdata.Protection Class (IEC 61140)</b>	1
<b>hei.techdata.Protection Class (IEC 60529)</b>	Device: IP 20 / Cable heating bath: IP 67